

Datasheet

Hepatitis C virus NS-3 monoclonal antibody, clone 12-15

Catalog Number: MAB6328

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against partial recombinant hepatitis C virus NS-3.

Clone Name: 12-15

Immunogen: Recombinant protein corresponding to least 60 amino acids in length of Chinese HCV NS-3.

Host: Mouse

Reactivity: Viruses

Applications: ELISA, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody reacts with synthetic HCV NS-3 Protein (amino acids 1378-1458).

Form: Lyophilized

Isotype: IgG1, kappa

Recommend Usage: Western Blot (1 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from 10 mM PBS, pH 7.0

Storage Instruction: Store at -20°C.

Aliquot after reconstitution to avoid repeated freezing and thawing.