

Datasheet

Hepatitis C virus NS-3 monoclonal antibody, clone 1-5

Catalog Number: MAB6772

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against partial recombinant Hepatitis C virus NS-3.

Clone Name: 1-5

Immunogen: Recombinant protein corresponding to least 60 residues in length of Chinese HCV NS-3.

Host: Mouse

Reactivity: Human

Applications: ELISA, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody has no cross reaction with HCV capsid region and other non-structural regions can be seen.

Form: Lyophilized

Isotype: IgG1, kappa

Recommend Usage: Western Blot (4 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from 10 mM PBS, pH 7.2

Storage Instruction: Store at -20°C.

Aliquot after reconstitution to avoid repeated freezing and thawing.