

## Datasheet

### **BPH monoclonal antibody, clone YPBH-2**

**Catalog Number:** MAB6778

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against Benign Prostatic Hyperplasia (BPH).

**Clone Name:** YPBH-2

**Immunogen:** Highly purified Benign Prostatic Hyperplasia (BPH).

**Host:** Mouse

**Reactivity:** Human

**Applications:** ELISA, IHC

(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** This antibody reacts with BPH or human prostate epithelial cells. It was nor reaction in non-prostate normal tissue, neither prostate malignant tissue positive. This antibody can react with some of well-moderately differential prostate normal tissue , and very rarely react with poor differential prostate normal tissue.

**Form:** Lyophilized

**Isotype:** IgG, kappa

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** Lyophilized from 10 mM PBS, pH 7.2

**Storage Instruction:** Store at -20°C.

Aliquot after reconstitution to avoid repeated freezing and thawing.