

Datasheet

CX3CL1 polyclonal antibody

Catalog Number: PAB16191

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against recombinant CX3CL1.

Immunogen: Recombinant protein corresponding to human CX3CL1.

Host: Rabbit

Reactivity: Human

Applications: ELISA

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Lyophilized

Recommend Usage: ELISA (0.5-1 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from PBS.

Storage Instruction: Store at -20°C on dry atmosphere.

After reconstitution with deionized water, store at -20°C or lower.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 6376

Gene Symbol: CX3CL1

Gene Alias: ABCD-3, C3Xkine, CXC3, CXC3C, NTN, NTT, SCYD1, fractalkine, neurotactin