

Datasheet

Beta-Galactosidase monoclonal antibody, clone DC1 4C7

Catalog Number: MAB7238

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against full length Beta-Galactosidase.

Clone Name: DC1 4C7

Immunogen: Native from bacteria Beta-Galactosidase.

Host: Mouse

Reactivity: Bacteria

Applications: ELISA, ICC, IF, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Beta-galactosidase (DC1 4C7).

Form: Liquid

Isotype: IgG1

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS

Storage Instruction: Store at 4°C for short term. For long term storage store at -20°C.
Aliquot to avoid repeated freezing and thawing.

References:

1. A sensitive micro-immunoassay using beta-galactosidase/anti-beta-galactosidase complexes. Durbin H, Bodmer WF. J Immunol Methods. 1987 Feb 26;97(1):19-27.