

Datasheet

Cdk2 monoclonal antibody, clone C15.2

Catalog Number: MAB7260

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against full length Cdk2.

Clone Name: C15.2

Immunogen: Full length xenopus Cdk2.

Host: Mouse

Reactivity: Xenopus

Applications: WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Cdk2 (C15.2). Predicted molecular weight 34 KDa.

Form: Liquid

Isotype: IgG2a

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS

Storage Instruction: Store at 4°C for short term. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.