

Datasheet

MST1R polyclonal antibody

Catalog Number: PAB16429

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of MST1R.

Immunogen: A synthetic peptide (conjugated with KLH) corresponding to human MST1R.

Host: Rabbit

Reactivity: Human

Applications: IHC-P

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Internal domain of human.

Form: Liquid

Purification: Immunoaffinity purification

Recommend Usage: Immunohistochemistry

(Formalin/PFA-fixed paraffin-embedded sections) (3 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.1% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -80°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 4486

Gene Symbol: MST1R

Gene Alias: CD136, CDw136, PTK8, RON