

Datasheet

PCGF1 polyclonal antibody (A01)

Catalog Number: H00084759-A01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant PCGF1.

Immunogen: PCGF1 (NP_116062, 105 a.a. ~ 187 a.a) partial recombinant protein with GST tag.

Sequence:

DSEEKRIREFYQSRGLDRVTQPTGEEPALSNLGLPFS
SFDHSKAHYYRYDEQLNLCLERLSSGKDKNKSVLQNK
YVRCSVRAE

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Ce, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 84759

Gene Symbol: PCGF1

Gene Alias: 2010002K04Rik, FLJ43754, MGC10882, NSPC1, RNF3A-2, RNF68

Gene Summary: PCGF1 is a mammalian homolog of the Drosophila polycomb group genes, which act as transcriptional repressors to regulate anterior-posterior patterning in early embryonic development (Nunes et al., 2001 [PubMed 11287196]). See also PCGF2 (MIM 600346).[supplied by OMIM]