

Datasheet

ZNF577 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00084765-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human ZNF577 protein.

Immunogen: ZNF577 (NP_116068, 1 a.a. ~ 478 a.a)
full-length human protein.

Sequence:

MSVRREQGSSSGEGSLSFEDVAVGFTREEWQFLDQS
QKVLYKEVMLENYINLVSIGYRGTKPDSLFLKEQGEPP
GIVEGAAHSQICPGFVIQSRRYAGKSDAFGGYGRSC
LHIKCDKLTGTGVKYHRCVKPSSPKSQLNDLQKICAGGK
PHECSVCGRAFSRKAQLIQHQRTERGEKPHGCGECG
KTFMRKIQLTEHQRTHTEGKPEHCSECGKAFSRKSQ
MVHQRTHTEGKPYRCSECGKAFSRKCRNLNRHQRSHT
GEKLYGCSVCGKAFSQKAYLTAHQRLHTGDKPYKCS
DCGRTFYFKSDLTRHQRIHTGEKPYECSECEKAFRSK
SKLIQHQRTHTEGERPYSCRECGKAFAHMSVLIKHEKT
HIRETAINSLTVEKPSSRSHTSLYMSELIQEQTNTVP
IEMPSSGTPPLLKSERLVGRNVVIVEQPFPRNQAFVV
NQEFEQRISLTNEVNVAPSVINYILYLTDIVSE

Host: Mouse

Reactivity: Human

Applications: WB-Ti, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 84765

Gene Symbol: ZNF577

Gene Alias: FLJ32291, MGC4400