

Datasheet

NFATC2IP MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00084901-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human NFATC2IP protein.

Immunogen: NFATC2IP (AAH21551, 1 a.a. ~ 138 a.a) full-length human protein.

Sequence:

MSEPLQSVVDHMAHLGVSPSRILLFGETELSPTATP
RTLKLGVADIIDCVVLTSSPEATETSQQLQLRVQGKEK
HQTLEVSLSRDSPLKTLMSHYEEAMGLSGRKLSFFFD
GTKLSGRELPADLGMESGDLIEVWG

Host: Mouse

Reactivity: Human

Applications: IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 84901

Gene Symbol: NFATC2IP

Gene Alias: FLJ14639, MGC126790, MGC138387