

Datasheet

TNS4 monoclonal antibody (M02), clone 1C1

Entrez GeneID: 84951**Gene Symbol:** TNS4**Gene Alias:** CTEN, FLJ14950, PP14434**Catalog Number:** H00084951-M02**Regulation Status:** For research use only (RUO)**Product Description:** Mouse monoclonal antibody
raised against a full-length recombinant TNS4.**Clone Name:** 1C1**Immunogen:** TNS4 (AAH18706, 1 a.a. ~ 480 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.**Sequence:**

MGSKASSPHGLGSPLVASPRLEKRLGGLAPQRGSRIS
VLSASPVSDVSYMFGSSQSLHSSNSSHQSSSRLES
PANSSSSLSLGSVSLCTRPSDFQAPRNPTLTMGQPR
TPHSPPLAKEHASSCPPSITNSMVDIPIVLINGCPEPGS
SPPQRTPGHQNSVQPGAASPSNPCPATRSNSQTLSD
APFTTCPEGPARDMQPTMKFVMDTSKYWFKPNITRE
QAIELLRKEEPGAFVIRDSSSYRGSFGLALKVQEVPA
AQSRPGEDSNDLIRHFLIESSAKGVHLKGADDEEYFGS
LSAFVCQHSIMALALPCKLTIPQRELGGADGASDSTDS
PASCQKKSAGCHTLYLSSVSVETLTGALAVQKAISTTF
ERDILPTPTVVHFKVTEQGITLTDVQRKVFFRRHYPLTT
LRFCGMDPEQRKWQKYCKPSWIFGFVAKSQTEPQEN
VCHLFAEYDMVQPASQVIGLVTALLQDAERM

Host: Mouse**Reactivity:** Human**Applications:** ELISA, WB-Re, WB-Tr(See our web site product page for detailed applications
information)**Protocols:** See our web site at<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols**Isotype:** IgG2a Kappa**Storage Buffer:** In 1x PBS, pH 7.4**Storage Instruction:** Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.