

Datasheet

LGALS12 purified MaxPab mouse polyclonal antibody (B02P)

Catalog Number: H00085329-B02P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human LGALS12 protein.

Immunogen: LGALS12 (NP_149092.2, 1 a.a. ~ 336 a.a) full-length human protein.

Sequence:

MSQPSGGRAPGTRIYSWSCPTVMSPGEKLDPIPDSFI
LQPPVFHPVVPYVTTIFGGLHAGKMVMLQGVPPLDAH
RFQVDFQCGCSLCPRPDIAFHFNPRFHTTKPHVICNTL
HGGRWQREARWPHLALRRGSSFLILFLFGNEEVKVS
NGQHFLHFRYRLPLSHVDTLGIFGDILVEAVGFLNINPF
VEGSREYPAGHPFLLMSPRLEVPCSHALPQGLSPGQV
IIVRGLVLQEPKHFTVSLRDQAAHAPVTLRASFADRTLA
WISRWGQKKLISAPFLFYPQRRFEVLLLFQEGGLKLAL
NGQGLGATSMNQQALEQLRELRISSVQLYCVHS

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 85329

Gene Symbol: LGALS12

Gene Alias: GALECTIN-12, GRIP1