

Datasheet

IgG (kappa) monoclonal antibody, clone KP-53

Catalog Number: MAB7930

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against native IgG, kappa.

Clone Name: KP-53

Immunogen: Native purified IgG, kappa light chain from
Bence Jones' gammopathy.

Host: Mouse

Reactivity: Human

Applications: IHC-P

(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Form: Liquid

Isotype: IgG2b

Recommend Usage: The optimal working dilution
should be determined by the end user.

Storage Buffer: In PBS, pH 7.3 (0.09% sodium azide)

Storage Instruction: Store at 4°C.