

Datasheet

Pentosidine monoclonal antibody, clone PEN-12

Catalog Number: MAB7939

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against synthetic peptide of Pentosidine.

Clone Name: PEN-12

Immunogen: A peptide (conjugated with HAS)
corresponding to Pentosidine .

Host: Mouse

Reactivity: Bovine,Dog,Human,Pig,Rat

Applications: ELISA, IHC

(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Form: Liquid

Concentration: 0.25 mg/mL

Isotype: IgG1

Recommend Usage: ELISA (0.1-0.5 ug/mL)

Immunohistochemistry (2 ug/mL)

The optimal working dilution should be determined by
the end user.

Storage Buffer: In ascites (50% Block Ace, 0.1%
proclin)

Storage Instruction: Store at -20°C.

Aliquot to avoid repeated freezing and thawing.