

Datasheet

Infectious bursal disease virus monoclonal antibody, clone Infectious bursal disease virus99

Catalog Number: MAB8522

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against Infectious bursal disease virus (IBDV).

Clone Name: IBDV99

Immunogen: Native purified very virulent IBD virus (field strain) derived from infected chicken bursas.

Host: Mouse

Reactivity: Chicken

Applications: ELISA, IHC, S-ELISA, WB
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: VP3 protein of IBD virus. MAbs detect IBD virus with high specificity and do not interact with other avian viruses.

Form: Liquid

Isotype: IgG1

Recommend Usage: Recommended pairs for IBD virus in sandwich immunoassay (capture-detection) :
IBDV99-IBDV105.

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (0.1% sodium azide)

Storage Instruction: Store at 4°C. Do not freeze.