

Datasheet

Influenza A virus Haemagglutinin H7 monoclonal antibody, clone 6B5

Catalog Number: MAB8531

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against Influenza A haemagglutinin H7 (H7N1).

Clone Name: 6B5

Immunogen: Influenza virus A/FPV (H7N1).

Host: Mouse

Reactivity: Viruses

Applications: HI

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Specificity: MAb 6B5 has low cross-reactivity to H1.

Form: Liquid

Isotype: IgG2b

Recommend Usage: The optimal working dilution
should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (0.1% sodium azide)

Storage Instruction: Store at 4°C. Do not freeze.