

## Datasheet

### Influenza A virus monoclonal antibody, clone 9F10

**Catalog Number:** MAB8541

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against Influenza A virus.

**Clone Name:** 9F10

**Immunogen:** Native nonstructural (NS) protein of influenza A virus.

**Host:** Mouse

**Reactivity:** Viruses

**Applications:** ELISA

(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** Non structural protein of influenza A virus.

**Form:** Liquid

**Isotype:** IgG2b

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (0.1% sodium azide)

**Storage Instruction:** Store at 4°C. Do not freeze.