

Datasheet

KLH monoclonal antibody, clone 11E8F5C3

Catalog Number: MAB8544

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against KLH .

Clone Name: 11E8F5C3

Immunogen: Keyhole Limpet Hemocyanin, KLH.

Host: Mouse

Applications: ELISA, WB

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Specificity: This antibody detected specifically
individual KLH and conjugated protein, as well as
KLH-conjugated peptide in both ELISA and Western
Blot.

Form: Lyophilized

Isotype: IgM

Recommend Usage: ELISA (1:2000-1:5000)

Western Blot (2-5 ug/mL)

The optimal working dilution should be determined by
the end user.

Storage Buffer: Lyophilized from PBS, pH 7.5, 0.1%
BSA.

Storage Instruction: Store at -20°C.

Aliquot after reconstitution to avoid repeated freezing
and thawing.