

Datasheet

SLIC1 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00124460-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human SLIC1 protein.

Immunogen: SLIC1 (ENSP00000300590, 1 a.a. ~ 129 a.a) full-length human protein.

Sequence:

MASPEHPGSPGCMGPITQCTARTQQEAPATGPDLP
PGPDGHLDTHTSGLSSNSSMTTRELQYWNQKCRW
KHVKLLFEIASARIEERKVSFVMGKSRPGEMTPGSR
GETGTAPEPDPRCPRQSDML

Host: Mouse

Reactivity: Human

Applications: Det Ab, IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 124460

Gene Symbol: SNX20

Gene Alias: MGC35578, SLIC-1, SLIC1