

Datasheet

MSI2 monoclonal antibody (M10), clone 3F11

Catalog Number: H00124540-M10

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a full length recombinant MSI2.

Clone Name: 3F11

Immunogen: MSI2 (AAH01526, 1 a.a. ~ 328 a.a)
full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

MEANGSQGTSGSANDSQHDPGKMFIGGLSWQTSPD
SLRDYFSKFGEIRECMVMRDPTTKRSRGFGVTFADP
ASVDKVLGQPHHELDSTIDPKVAFPRRAQPKMVRT
KKIFVGGLSANTVVEDVKQYFEQFGKVEDAMLMFDKT
TNRHRGFGFVTFENEDVVEKVCEIHFHEINNKMVECK
KAQPKVEMFPPGTRGRARGLPYTMDAFMLGMGMLG
YPNFVATYGRGYPGFAPSYGYQFPGFPAAAYGPVAA
AAVAAARGSGSNPARPGGFPGANSPGPVADLYGPAS
QDSGVGNYISAASPPGSGFGHGIAGPLIATAFTNGY
H

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA, WB-Re
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 124540

Gene Symbol: MSI2

Gene Alias: FLJ36569, MGC3245, MSI2H

Gene Summary: This gene encodes a protein containing two conserved tandem RNA recognition motifs. Similar proteins in other species function as RNA-binding proteins and play central roles in posttranscriptional gene regulation. Two transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq]