

Datasheet

CANT1 MaxPab rabbit polyclonal antibody (D01)

Catalog Number: H00124583-D01

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human CANT1 protein.

Immunogen: CANT1 (NP_620148.1, 1 a.a. ~ 401 a.a) full-length human protein.

Sequence:

MPVQLSEHPEWNE SMHSLRISVGGLPVLASMTKAADP
RFRPRWKVILTFVGAAILWLLCSHRPAPGRPPTHNAH
NWRLGQAPANWYNDTYPLSPQRTYPAGIRYRIAVIAD
LDTESRAQEENTWFSYLKKGYLTLSDSGDKVAVEWD
KDHGVLESHLAEKGRGMELSDLIVFNGKLYSVDDRTG
VVYQIEGSKAVPWVILSDGDGTVEKGFKAEWLAVKDE
RLYVGGLGKEWTTTTGDVVNENPEWVKVVGKGSVD
HENWVSNYNALRAAAGIQPPGYLIHESACWSDTLQR
WFFLPRRASQERYSEKDDERKGANLLLSASPDFGDIA
VSHVGAVVPTHGFSSFKFIPNTDDQIIVALKSEEDSGR
VASYIMAFITLDGRFLLPETKIGSVKYEIEFI

Host: Rabbit

Reactivity: Human, Mouse

Applications: IP, WB-Ce, WB-Ti, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 124583

Gene Symbol: CANT1

Gene Alias: SCAN-1, SHAPY