

Datasheet

ZBPB2 MaxPab mouse polyclonal antibody (B01)

interaction. Dun MD, Smith ND, Baker MA, Lin M, Aitken RJ, Nixon B. J Biol Chem. 2011 Aug 31. [Epub ahead of print]

Catalog Number: H00124626-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human ZBPB2 protein.

Immunogen: ZBPB2 (ENSP00000335384, 1 a.a. ~ 315 a.a) full-length human protein.

Sequence:

MRTCVLLSAVLWCLTGDKIYVELHQNSPVLICMDFKLS
KKEIVDPTYLWIGPNEKTLTGNNRINITETGQLMVKDFL
EPLSGLYTCTLSYKTVKAETQEEKTVKKRYDFMVFAY
REPDSYQMAVRFTTRSCIGRYNDVFFRVLKKILDSLIS
DLSCHVIEPSYKCHSVEIPEHGLIHELFIQVNPFPAPG
WKGACNGSVDCEDTTNHNILQARDRIEDFFRSQAYIF
YHNFNKTLPAMHFVDHSLQVVRLDSCRPGFGKNERL
HSNCASCCVVCSPATFSPDVNVTCQTCVSVLTYGAKS
CPQTSNKNQQYED

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 124626

Gene Symbol: ZBPB2

Gene Alias: MGC41930, ZPBPL

References:

1. The chaperonin containing TCP1 complex (CCT/TRiC) is involved in mediating sperm-oocyte