

Datasheet

TRAPPC5 monoclonal antibody (M01), clone 3D9

Catalog Number: H00126003-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a full length recombinant TRAPPC5.

Clone Name: 3D9

Immunogen: TRAPPC5 (AAH42161, 1 a.a. ~ 188 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

MEARFTRGKSALLERALARPRTEVSLSAFALLFSELVQ
HCQSRVFSVAELQSRLAALGRQVGARVLDALVAREKG
ARRETKVLGALLFVKGAVWKALFGKEADKLEQANDDA
RTFYIIEREPLINTYISVPKENSTLNCASFTAGIVEAVLTH
SGFPAKVTAHWHKGTTLMIKFEEAVIARDRALEGR

Host: Mouse

Reactivity: Human

Applications: ELISA

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 126003

Gene Symbol: TRAPPC5

Gene Alias: MGC52424