

Datasheet

NTN5 purified MaxPab mouse polyclonal antibody (B01P)

Gene Alias: -

Catalog Number: H00126147-B01P

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human NTN5 protein.

Immunogen: NTN5 (NP_665806.1, 1 a.a. ~ 489 a.a) full-length human protein.

Sequence:

MPVTFALLLLGQATADPCYDPQGRPQFCLPPVTQLA
AVAASCPQACALSPGNHLGARETCNGSLTLALGGPFL
LTSVSLRFCTPGPPALILSAAWASGGPWRLWHRPAW
PGALGGPERVTFHSTPGPKATVAASHLRVEFGGQAGL
AAAGLRGRCQCHGHAARCAARARPPRCHCRHHTTG
PGCESCRPSHRDWPWRPATPRHPCLPCSCNQHA
RRCRFNSELFRLSGGRSGGVCERCRRHHTAGRHCHYC
QPGFWRDPSQPIFSRRACRACQCHPIGATGGTCNQ
SGQCTCKLGVTGLTCNRCPGYQQSRSPRMPQCQRIP
EATTTLATTPGAYSSDPQCQNYCNMSDTRVHMSLRR
YCQQDHVLAQVLASEAAGPAWQRLAVRVLAVYKQR
AQPVRRGDQDAWVPRADLTGCLRLQPGTDYLLLS
AVGDPDPTRLILDRHGLALPWRPRWARPLKRLQQEER
AGGCRGVRAPTPSPRPEH

Host: Mouse

Reactivity: Human

Applications: IF, WB-Ce, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 126147

Gene Symbol: NTN5