

Datasheet

MGC29891 monoclonal antibody (M01), clone 3A3-1F4

Gene Symbol: GABPB2

Gene Alias: MGC29891

Catalog Number: H00126626-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full length recombinant MGC29891.

Clone Name: 3A3-1F4

Immunogen: MGC29891 (AAH27033, 1 a.a. ~ 448 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MSLVDLGKRLLEAARKGQDDEVRTL MANGAPFTTDW
LGTSPLHLAAQYGHYSTAEVLLRAGVSRDARTKVDRT
PLHMAAADGHAHIVELLVRNGADVNAK DMLKMTALH
WATERHHRDVVELLIKYGADVHAFSKFDKSAFDIALEK
NNAEILVILQEAMQNQVNVNPERANPVTDPVSMAAPFI
FTSGEVVNLASLISSTNTKTTSGDPHASTVQFSNSTTS
VLATLAALAEASVPLSNSHRATANTEEIIIEGNSVDSSIQ
QVMGSGGQRVITIVTDGVPLGNIQTSIPTGGIGQPFIVT
VQDGQQVLTVPAGKVAEETVIKEEEEEKLPLTKKPRIG
EKTNSVEESKEGNERELLQQQLQEANRR AQEYRHQL
LKKEQEAEQYRLKLEAIARQQPNGVDFTMV EEEVAEVD
AVVVTEGELEERETKVTGSAGTTEPHTRVSMATVSS

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA, WB-Re, WB-Tr
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2b kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 126626