

Datasheet

TMEM42 MaxPab mouse polyclonal antibody (B01)

Catalog Number: H00131616-B01

Regulation Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human TMEM42 protein.

Immunogen: TMEM42 (NP_653239, 1 a.a. ~ 159 a.a) full-length human protein.

Sequence:

MAERPGPPGGAVSATAYPDTPAEFPFHLQAGAMRRR
FWGVFNCLCAGAFGALAAASAKLAFGSEVSMGLCVL
GIIVMASTNSLMWTFFSRGLSFSMSSAIASVTVTFSNIL
SSAFLGYVLYGECQEVLLWWGGVFLILCGLTLIHRKLPP
TWKPLPHKQQ

Host: Mouse

Reactivity: Human

Applications: Det Ab, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: No additive

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 131616

Gene Symbol: TMEM42

Gene Alias: MGC29956