

Datasheet

POU5F2 monoclonal antibody (M03), clone 1B12

Gene Alias: DKFZp686P02123, FLJ25680, SPRM-1

Catalog Number: H00134187-M03

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full-length recombinant POU5F2.

Clone Name: 1B12

Immunogen: POU5F2 (AAH29532, 1 a.a. ~ 304 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MPLRVDTLTWLSTQAAPGRVMVWPAVRPGICPGPDV
WRIPLGPLPHEFRGWIAPCRPRLGASEAGDWLRRPSE
GALPGPYIALRSIPKLPPPEDISGILKELQQLAKELRQK
RLSLGYSQADVGIavgALFGKVLSTTICRFEAQQLSV
ANMWKLRPLLKKWLKEVEAENLLGLCKMEMILQQSGK
WRRASRERRIGNSLEKFFQRCPKPTPQQISHIAGCLQL
QKDVVRVWFYNRSKMGSRPTNDASPREIVGTAGPPC
PGAPVCFHLGLGLPVDIPHYTRYLYSAGVAHSSAPATTL
GLLRF

Host: Mouse

Reactivity: Human

Applications: ELISA

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG1 Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 134187

Gene Symbol: POU5F2