

Datasheet

IgG (Fc) monoclonal antibody, clone EM-07 (FITC)

Catalog Number: MAB9827

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against IgG (Fc).

Clone Name: EM-07

Immunogen: Fusion protein of human IgG Fc fragment.

Host: Mouse

Reactivity: Human

Applications: IF

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Form: Liquid

Conjugation: FITC

Concentration: 0.1 mg/mL

Isotype: IgG1

Recommend Usage: Immunofluorescence (2-5 ug/mL)

The optimal working dilution should be determined by
the end user.

Storage Buffer: In PBS, pH 7.4 (0.02% BSA, 15 mM
sodium azide)

Storage Instruction: Store in the dark at 4°C. Avoid
prolonged exposure to light. Do not freeze.