

Datasheet

Itgb2 monoclonal antibody, clone C71/16 (Biotin)

Catalog Number: MAB9870

Regulation Status: For research use only (RUO)

Product Description: Rat monoclonal antibody raised against native Itgb2.

Clone Name: C71/16

Immunogen: Native Itgb2 from cell membrane of mouse T cell lymphoma BW5147.

Host: Rat

Reactivity: Mouse

Applications: Flow Cyt

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Conjugation: Biotin

Isotype: IgG2a, kappa

Recommend Usage: Flow Cytometry (? 1 ug/10⁶ cells)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.09% sodium azide)

Storage Instruction: Store at 4°C. Do not freeze.

Entrez GeneID: 16414

Gene Symbol: Itgb2

Gene Alias: 2E6, AI528527, Cd18, LAD, LCAMB, Lfa1, MF17