

## Datasheet

### FAM59B polyclonal antibody

**Catalog Number:** PAB16861

**Regulation Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of FAM59B.

**Immunogen:** A synthetic peptide corresponding to C-terminus 19 amino acids of human FAM59B.

**Host:** Rabbit

**Reactivity:** Human, Mouse, Rat

**Applications:** ELISA, IHC, WB-Ce  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Recommend Usage:** Western Blot (1-2 ug/mL)  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS (0.02% sodium azide)

**Storage Instruction:** Store at 4°C for three months. For long term storage store at -20°C.  
Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 150946

**Gene Symbol:** FAM59B

**Gene Alias:** FLJ00375, KIAA2038