

Datasheet

Kcnp1 polyclonal antibody

Catalog Number: PAB16878

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of Kcnp1.

Immunogen: A synthetic peptide (conjugated with KLH) corresponding to mouse Kcnp1.

Sequence: QTKQRRPSKDK

Host: Rabbit

Reactivity: Mouse

Applications: EIA, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Lyophilized

Recommend Usage: Enzyme Immunoassay (1:300000)

Western Blot (1:4000)

The optimal working dilution should be determined by the end user.

Storage Buffer: Lyophilized from 10 mM PBS

Storage Instruction: Store at -20°C on dry atmosphere.

After reconstitution with deionized water, store at -20°C or lower.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 70357

Gene Symbol: Kcnp1

Gene Alias: 2900046L02Rik, 3202002F18Rik, A1842357, AW494471, KCHIP1, Kchip1.2