

Datasheet

C-peptide monoclonal antibody, clone 7E10

Catalog Number: MAB2807

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against synthetic peptide of C-peptide.

Clone Name: 7E10

Immunogen: A synthetic peptide corresponding to
human C-peptide.

Host: Mouse

Reactivity: Human

Applications: C-ELISA, EIA, RIA

(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Specificity: This antibody bind with equal efficiency to
the human free C-peptide and C-peptide region in
proinsulin molecules.

Form: Liquid

Isotype: IgG1, kappa

Recommend Usage: The optimal working dilution
should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (0.09% sodium azide)

Storage Instruction: Store at 4°C.