

Datasheet

Ghrelin monoclonal antibody, clone Ghr5

Catalog Number: MAB2827

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against Ghrelin.

Clone Name: Ghr5

Immunogen: Ghrelin.

Host: Mouse

Applications: ELISA

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Isotype: IgG2a

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (0.09% sodium azide)

Storage Instruction: Store at 4°C.