

Datasheet

Human chorionic gonadotropin monoclonal antibody, clone F1

Catalog Number: MAB2830

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against human chorionic gonadotropin (hCG).

Clone Name: F1

Immunogen: HCG.

Host: Mouse

Reactivity: Human

Applications: ELISA

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Specificity: This antibody is specific for the alpha-chain
of FSH, LH, TSH and HCG.

Form: Liquid

Isotype: IgG

Recommend Usage: The optimal working dilution
should be determined by the end user.

Storage Buffer: In 10 mM Tris, 150 mM NaCl, pH 7.5
(0.05% sodium azide)

Storage Instruction: Store at 4°C.