

Datasheet

Human chorionic gonadotropin monoclonal antibody, clone 76G3

Catalog Number: MAB2836

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against human chorionic gonadotropin (hCG).

Clone Name: 76G3

Immunogen: HCG.

Host: Mouse

Reactivity: Human

Applications: S-ELISA

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Specificity: This antibody specifically recognizes HCG;
it reacts mainly with alpha-subunit. Cross-reacting with
whole molecule of HCG, LH, TSH and FSH.
Cross-reactivity with free beta-subunit of HCG was not
detected.

Form: Liquid

Isotype: IgG2b

Recommend Usage: The optimal working dilution
should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (0.09% sodium azide)

Storage Instruction: Store at 4°C.