

Datasheet

Insulin/Proinsulin monoclonal antibody, clone D6C4

Catalog Number: MAB2844

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against Insulin/Proinsulin.

Clone Name: D6C4

Immunogen: Rat and mouse Insulin/Proinsulin.

Host: Mouse

Reactivity: Bovine, Human, Mouse, Pig, Rat

Applications: EIA, IHC-Fr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody cross react with insulin and proinsulin.

Form: Liquid

Isotype: IgG1, kappa

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (0.09% sodium azide)

Storage Instruction: Store at 4°C.