

Datasheet

Progesterone monoclonal antibody, clone XM207

Catalog Number: MAB2858

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against Progesterone.

Clone Name: XM207

Immunogen: Progesterone conjugated with BSA.

Host: Mouse

Applications: ELISA

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Specificity: There is no cross-reactivity with BSA.

Form: Liquid

Isotype: IgG2b

Recommend Usage: The optimal working dilution
should be determined by the end user.

Storage Buffer: In 10 mM Tris, 150 mM NaCl, pH 7.5
(0.05% sodium azide)

Storage Instruction: Store at 4°C.