

Datasheet

IFNA2 monoclonal antibody, clone C10F5

Catalog Number: MAB3001

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against recombinant IFNA2.

Clone Name: C10F5

Immunogen: Recombinant protein corresponding to
human IFNA2.

Host: Mouse

Reactivity: Human

Applications: S-ELISA

(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Form: Liquid

Isotype: IgG1

Recommend Usage: The optimal working dilution
should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (0.09% sodium azide)

Storage Instruction: Store at 4°C.

Entrez GeneID: 3440

Gene Symbol: IFNA2

Gene Alias: IFNA, INFA2, MGC125764, MGC125765