

## Datasheet

### Procalcitonin monoclonal antibody, clone 27A3

**Catalog Number:** MAB3017

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody  
raised against partial Procalcitonin.

**Clone Name:** 27A3

**Immunogen:** Human Procalcitonin conjugated with  
carrier protein.

**Host:** Mouse

**Reactivity:** Human

**Applications:** ELISA, WB

(See our web site product page for detailed applications  
information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product  
page for detailed protocols

**Form:** Liquid

**Isotype:** IgG2a

**Recommend Usage:** The optimal working dilution  
should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (0.09% sodium azide)

**Storage Instruction:** Store at 4°C.