

## Datasheet

### **S100A1/S100B monoclonal antibody, clone 6G1**

**Catalog Number:** MAB3029

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against native S100A1/S100B.

**Clone Name:** 6G1

**Immunogen:** Native purified human brain S100A1/S100B.

**Host:** Mouse

**Reactivity:** Human

**Applications:** S-ELISA, WB

(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** This antibody is specific to S100 beta-beta and S100 alpha-beta.

**Form:** Liquid

**Isotype:** IgG1

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS, pH 7.4 (0.09% sodium azide)

**Storage Instruction:** Store at 4°C.

**Entrez GeneID:** 6271|6285

**Gene Symbol:** S100A1

**Gene Alias:** S100, S100-alpha, S100A