

Datasheet

S100A1/S100B monoclonal antibody, clone 3B10

Catalog Number: MAB3030

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against native S100A1/S100B.

Clone Name: 3B10

Immunogen: Native purified human brain S100A1/S100B.

Host: Mouse

Reactivity: Human

Applications: S-ELISA, WB

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody is specific to S-100 alpha-beta.

Form: Liquid

Isotype: IgG2a

Recommend Usage: The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, pH 7.4 (0.09% sodium azide)

Storage Instruction: Store at 4°C.

Entrez GeneID: 6271|6285

Gene Symbol: S100A1

Gene Alias: S100, S100-alpha, S100A