

Datasheet

Tlr11 polyclonal antibody

Catalog Number: PAB17961

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of Tlr11.

Immunogen: A synthetic peptide corresponding to the residues surrounding amino acid 793 of mouse Tlr11.

Host: Rabbit

Reactivity: Human, Mouse, Rat

Applications: WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody recognizes ~70-90 KDa of TLR11 in Jurkat cell lysate.

Form: Liquid

Recommend Usage: Western Blot (0.5-4 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS pH 7.2 (30% glycerol, 0.5% BSA, 0.01% thimerosal)

Storage Instruction: Store at -20°C. For long term storage store at -80°.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 239081

Gene Symbol: Tlr11

Gene Alias: Gm287, MGC129352, MGC129353