

Datasheet

ITPR1 polyclonal antibody

Gene Symbol: ITPR1

Catalog Number: PAB18402

Gene Alias: INSP3R1, IP3R, IP3R1, SCA15, SCA16

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of ITPR1.

Immunogen: A synthetic peptide corresponding to human ITPR1.

Sequence: A synthetic peptide corresponding to amino acids 1565-1615 of human ITPR1.

Host: Rabbit

Reactivity: Human, Mouse, Rat

Applications: ELISA, IF, IHC-P
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: This antibody is specific to ITPR1.

Form: Liquid

Purification: Affinity purification

Concentration: 1 mg/mL

Recommend Usage: Immunohistochemistry
(1:50-1:100)
Immunofluorescence (1:500-1:1000)
ELISA (1:20000)
The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS, 150mM NaCl, pH 7.4 (50% glycerol, 0.02% sodium azide)

Storage Instruction: Store at -20°C.
Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 3708