

Datasheet

LIMD1 monoclonal antibody (M01), clone 2G5

Catalog Number: H00008994-M01

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant LIMD1.

Clone Name: 2G5

Immunogen: LIMD1 (NP_055055.1, 577 a.a. ~ 674 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

VDSENKIYCVRDYHKVLAPKCAACGLPILPPEGSDETI
RVVSMDRDYHVECYHCEDCGLELNDEDGHRCPLED
HLFCHSCHVKRLEKRPSSTALHQH

Host: Mouse

Reactivity: Human

Applications: ELISA, IF, IHC-P, S-ELISA, WB-Re,
WB-Tr
(See our web site product page for detailed applications
information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2b Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 8994

Gene Symbol: LIMD1

Gene Alias: -