

Datasheet

MPZL1 purified MaxPab rabbit polyclonal antibody (D01P)

Catalog Number: H00009019-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human MPZL1 protein.

Immunogen: MPZL1 (NP_003944.1, 1 a.a. ~ 269 a.a) full-length human protein.

Sequence:

MAASAGAGAVIAAPDSRRWLWSVLAAALGLLTAGVSA
LEVYTPKEIFVANGTQGKLTCKFKSTTTGGLTSVSWS
FQPEGADTTVSFFHYSQGQVYLGNYPPFKDRISWAG
DLDKKDASINIENMQFIHNGTYICDVKNPPDIVVQPGHI
RLYVVEKENLPVFPVWVVVGIVTAVVLGLTLLISMILAV
LYRRKNSKRDTGCGSTSESLSPVKQAPRKSPSDTEGL
VKSLPSGSHQGPVIYAQLDHSGGHHSDKINKSESVVY
ADIRKN

Host: Rabbit

Reactivity: Human

Applications: IF, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 9019

Gene Symbol: MPZL1

Gene Alias: FLJ21047, PZR, PZR1b, PZRa, PZRb