

Datasheet

SOCS3 monoclonal antibody (M03), clone 1F6

Catalog Number: H00009021-M03

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a full-length recombinant SOCS3.

Clone Name: 1F6

Immunogen: SOCS3 (AAH60858, 1 a.a. ~ 225 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

MVTHSKFPAAGMSRPLDTSLRLKTFSSKSEYQLVVNA
VRKLQESGFYWSAVTGGEANLLLSAEPAGTFLIRDSS
DQRHFFTL SVKTSQGTKNLRIQCEGGSFSLQSDPRST
QPVPRFDCVLKLVHHYMPPPGAPSFSPPTESSEVP
EQPSAQPLPGSPRRAYYIYSGGEKIPLVLSRPLSSNV
ATLQHLCRKTVNGHLDSEYKVTQLPGPIREFLDQYDA
PL

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 9021

Gene Symbol: SOCS3

Gene Alias: ATOD4, CIS3, Cish3, MGC71791,
SOCS-3, SSI-3, SSI3