

## Datasheet

### **BRSK2 purified MaxPab mouse polyclonal antibody (B01P)**

**Catalog Number:** H00009024-B01P

**Regulation Status:** For research use only (RUO)

**Product Description:** Mouse polyclonal antibody raised against a full-length human BRSK2 protein.

**Immunogen:** BRSK2 (AAH24291.1, 1 a.a. ~ 170 a.a) full-length human protein.

**Sequence:**

MSNLTPESSPELAKKSWFGNFISLEKEEQIFVVIKDKPL  
SSIKADIVHAFLSIPSLSHSVISQTSFRAEYKATGGPAVF  
QKPVKFQVDITYTEGGEAQKENGISVTFTLLSGPSRR  
FKRVVETIQAQLLSTHDPAAQHLSEPPPPAPGLSWG  
AGLKGQKVATSYESSL

**Host:** Mouse

**Reactivity:** Human

**Applications:** WB-Tr

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Storage Buffer:** In 1x PBS, pH 7.4

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 9024

**Gene Symbol:** BRSK2

**Gene Alias:** C11orf7, FLJ41362, PEN11B, SAD1, STK29