

Datasheet

HMG20A purified MaxPab rabbit polyclonal antibody (D01P)

Catalog Number: H00010363-D01P

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against a full-length human HMG20A protein.

Immunogen: HMG20A (NP_060670.1, 1 a.a. ~ 347 a.a) full-length human protein.

Sequence:

MENLMTSSTLPPLFADEDGSKESNDLATTGLNHPEVP
YSSGATSSTNNPEFVEDLSQGQLQSESSNAAEGNE
QRHEDEQRSKRGGWSKGRKRKKPLRDSNAPKSPLT
GYVRFMNERREQLRAKRPEVPFPEITRMLGNEWSKLP
PEEKQRYLDEADRDKERYMKELEQYQKTEAYKVFSR
KTQDRQKGKSHRQDAARQATHDHEKETEVKERSVFD
IPIFTEEFLNHSKAREALRQLRKSNMFEERNAAALQK
HVESMRTAVEKLEVDVIQERSRNTVLQQHLETLRQVL
TSSFASMPPLPGSGETPTVDTIDSYMNRLHSIILANPQD
NENFIATVREVVNRLDR

Host: Rabbit

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 10363

Gene Symbol: HMG20A

Gene Alias: FLJ10739, HMGX1, HMGXB1