

Datasheet

KLF2 monoclonal antibody (M08), clone 1D1

Catalog Number: H00010365-M08

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody
raised against a partial recombinant KLF2.

Clone Name: 1D1

Immunogen: KLF2 (NP_057354, 263 a.a. ~ 350 a.a)
partial recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.

Sequence:

SWPRKRTATHTCSYAGCGKTYTKSSHLKAHLRTHTGE
KPYHCNWDGCGWKFARSDDELTRHYRKHTGHRPFQC
HLCDFRFSRSDHLALH

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications
information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.

Entrez GeneID: 10365

Gene Symbol: KLF2

Gene Alias: LKLF