

Datasheet

TUBB4B monoclonal antibody (M02), clone 1G3

Entrez GeneID: 10383**Gene Symbol:** TUBB4B**Gene Alias:** RP13-122B23.2, Beta2, TUBB2, TUBB2C**Catalog Number:** H00010383-M02**Regulation Status:** For research use only (RUO)**Product Description:** Mouse monoclonal antibody
raised against a full length recombinant TUBB4B.**Clone Name:** 1G3**Immunogen:** TUBB4B (AAH01911, 1 a.a. ~ 445 a.a)
full-length recombinant protein with GST tag. MW of the
GST tag alone is 26 KDa.**Sequence:**

MREIVHLQAGQCGNQIGAKFWEVISDEHGIDPTGTYH
GDSDLQLERINVYYNEATGGKYVPRAVLVDLEPGTMD
SVRSGPFGQIFRPDNFVFGQSGAGNNWAKGHYTEGA
ELVDSVLDVVRKEAESCDCLQGFQLTHSLGGGTGSG
MGTLISKIREEYPDRIMNTFSVVPSPKVSDTVVEPYNA
TLSVHQLVENTDETYCIDNEALYDICFRTLKLTPTYGD
LNHLVSATMSGVTTCLRFPGLNADLRKLAVNMVPFP
RLHFFMPGFAPLTSRGSQQYRALTVPELTQQMFDAKN
MMAACDPRHGRYLTVAAVFRGRMSMKEVDEQMLNV
QKNSSYFVEWIPNNVKTAVCDIPPRGLKMSATFIGNS
TAIQELFKRISEQFTAMFRRKAFLHWYTGEGMDEMEF
TEAESNMNDLVSEYQQYQDATAEEEEGEFEEEEEEV
A

Host: Mouse**Reactivity:** Human**Applications:** ELISA, IF, S-ELISA, WB-Ce, WB-Re
(See our web site product page for detailed applications
information)**Protocols:** See our web site at
<http://www.abnova.com/support/protocols.asp> or product
page for detailed protocols**Isotype:** IgG1 Kappa**Storage Buffer:** In 1x PBS, pH 7.4**Storage Instruction:** Store at -20°C or lower. Aliquot to
avoid repeated freezing and thawing.